

Product Datasheet

Dynorphin A (2-13) (orb364576)

Description	Peptide to Dynorphin A (2-13) (Porcine)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Gly-Phe-Leu-Arg-Arg-Ile-Arg-Pro-Lys-Leu-Lys
Formula	$C_{66}H_{117}N_{23}O_{13}$
Purity	> 95%
Source	Synthetic
MW	1440.81
Expiration Date	12 months from date of receipt.

