

Product Datasheet

EAE CP (orb364587)

Description	Peptide to EAE CP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Val-Phe-Ile-Leu-Gly-Pro-Leu-Arg-Leu-Leu-Gly
Formula	$C_{59}H_{100}N_{14}O_{12}$
Purity	> 95%
Source	Synthetic
MW	1197.54
Expiration Date	12 months from date of receipt.

