

Product Datasheet

Ep-CAM (263-271) (orb364612)

Description	Peptide to Ep-CAM (263-271)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Leu-Lys-Ala-Gly-Val-Ile-Ala-Val
Formula	$C_{38}H_{70}N_{10}O_{10}$
Purity	> 95%
Source	Synthetic
MW	827.04
Expiration Date	12 months from date of receipt.