

Product Datasheet

EPE-Peptide (orb364613)

Description	Peptide to EPE-Peptide (Human)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Phe-Ser-Trp-Gly-Ala-Glu-Gly-Gln-Arg
Formula	$C_{46}H_{64}N_{14}O_{14}$
Purity	> 95%
Source	Synthetic
MW	1037.11
CAS Number	29705-92-8
Expiration Date	12 months from date of receipt.