

Product Datasheet

Fmoc-Lys(Boc)-Leu-Lys(Boc) (orb364634)

Description	Peptide to Fmoc-Lys(Boc)-Leu-Lys(Boc)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Fmoc-Lys (Boc)-Leu-Lys (Boc)
Formula	$C_{26}H_{35}N_8O_7S$
Purity	> 95%
Source	Synthetic
MW	809.90
Expiration Date	12 months from date of receipt.