

Product Datasheet

FMRF-related peptide, SDPFLRF-NH2 (orb364639)

Description	Peptide to FMRF-related peptide, SDPFLRF-NH2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Asp-Pro-Phe-Leu-Arg-Phe-NH2
Formula	$C_{42}H_{61}N_{11}O_{10}$
Purity	> 95%
Source	Synthetic
MW	880.02
Expiration Date	12 months from date of receipt.