

## Product Datasheet

### Forkhead derived peptide (orb364640)

|                         |                                                                         |
|-------------------------|-------------------------------------------------------------------------|
| <b>Description</b>      | Peptide to Forkhead derived peptide                                     |
| <b>Form/Appearance</b>  | Lyophilized powder                                                      |
| <b>Storage</b>          | Shipped at 4°C. Store at -20°C for one year.                            |
| <b>Note</b>             | For research use only                                                   |
| <b>Protein Sequence</b> | Lys-Lys-Ile-Ser-Gly-Arg-Leu-Ser-Pro-Ile-Met-Thr-Glu-Gln-NH <sub>2</sub> |
| <b>Formula</b>          | C <sub>68</sub> H <sub>123</sub> N <sub>21</sub> O <sub>20</sub> S      |
| <b>Purity</b>           | > 95%                                                                   |
| <b>Source</b>           | Synthetic                                                               |
| <b>MW</b>               | 1586.93                                                                 |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                         |