

Product Datasheet

Galanin (2-11) peptide (orb364647)

Description	Galanin (2-11) peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Trp-Thr-Leu-Asn-Ser-Ala-Gly-Tyr-Leu-Leu-NH ₂
Formula	C ₅₄ H ₈₁ N ₁₃ O ₁₄
Purity	> 95%
Source	Synthetic
MW	1136.33
CAS Number	367518-31-8
Expiration Date	12 months from date of receipt.

