

Product Datasheet

GIP (1-30) (orb364661)

Description	Peptide to GIP (1-30) (Porcine)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Ala-Glu-Gly-Thr-Phe-Ile-Ser-Asp-Tyr-Ser-Ile-Ala-Met-Asp-Lys-Ile-Arg-Gln-Gln-Asp-Phe-Val-Asn-Trp-Leu-Leu-Ala-Gln-Lys
Formula	$C_{162}H_{244}N_{40}O_{48}S$
Purity	> 95%
Source	Synthetic
MW	3552.10
Expiration Date	12 months from date of receipt.