

Product Datasheet

GIP (3-42) (orb364663)

Description	Peptide to GIP (3-42) (Human)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Glu-Gly-Thr-Phe-Ile-Ser-Asp-Tyr-Ser-Ile-Ala-Met- Asp-Lys-Ile-His-Gln-Gln-Asp-Phe-Val-Asn-Trp-Leu- Leu-Ala-Gln-Lys-Gly-Lys-Lys-Asn-Asp-Trp-Lys-His- Asn-Ile-Thr-Gln
Formula	C ₂₁₄ H ₃₂₄ N ₅₈ O ₆₃ S
Purity	> 95%
Source	Synthetic
MW	4749.38
Expiration Date	12 months from date of receipt.