

Product Datasheet

Glucagon peptide (orb364670)

Description	Glucagon peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	His-Ser-Gln-Gly-Thr-Phe-Thr-Ser-Asp-Tyr-Ser-Lys-Tyr-Leu-Asp-Ser-Arg-Arg-Ala-Gln-Asp-Phe-Val-Gln-Trp-Leu-Met-Asn-Thr-Lys-Arg-Asn-Lys-Asn-Asn-Ile-Ala
Formula	$C_{192}H_{295}N_{59}O_{60}S$
Purity	> 95%
Source	Synthetic
MW	4421.92
CAS Number	62340-29-8
Expiration Date	12 months from date of receipt.