

Product Datasheet

Glycogen Synthase-derived peptide (orb364677)

Description	Peptide to Glycogen Synthase-derived peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Lys-Leu-Asn-Arg-Thr-Leu-Ser-Val-Ala
Formula	$C_{49}H_{92}N_{16}O_{14}$
Purity	> 95%
Source	Synthetic
MW	1129.38
Expiration Date	12 months from date of receipt.