

Product Datasheet

GnT-V (nt38-67) (orb364687)

Description	Peptide to GnT-V (nt38-67)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Val-Leu-Pro-Asp-Val-Phe-Ile-Arg-Cys-Val
Formula	$C_{54}H_{89}N_{13}O_{13}S$
Purity	> 95%
Source	Synthetic
MW	1159.45
Expiration Date	12 months from date of receipt.