

Product Datasheet

gp100 (154-162) (orb364688)

Description	Peptide to gp100 (154-162)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Thr-Trp-Gly-Gln-Tyr-Trp-Gln-Val
Formula	$C_{58}H_{78}N_{14}O_{14}$
Purity	> 95%
Source	Synthetic
MW	1195.35
Expiration Date	12 months from date of receipt.

