

Product Datasheet

gp100 (177-186) (orb364690)

Description	Peptide to gp100 (177-186)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Met-Leu-Gly-Thr-His-Thr-Met-Glu-Val
Formula	$C_{45}H_{76}N_{12}O_{15}S_2$
Purity	> 95%
Source	Synthetic
MW	1089.31
Expiration Date	12 months from date of receipt.