

Product Datasheet

gp100 (178-187) (orb364691)

Description	Peptide to gp100 (178-187)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Met-Leu-Gly-Thr-His-Thr-Met-Glu-Val
Formula	$C_{42}H_{71}N_{11}O_{14}S_2$
Purity	> 95%
Source	Synthetic
MW	1018.23
Expiration Date	12 months from date of receipt.