

Product Datasheet

gp100 (209-217) (orb364692)

Description	Peptide to gp100 (209-217)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ile-Met-Asp-Gln-Val-Pro-Phe-Ser-Val
Formula	$C_{47}H_{74}N_{10}O_{14}S$
Purity	> 95%
Source	Synthetic
MW	1035.24
Expiration Date	12 months from date of receipt.

