

Product Datasheet

gp100 (476-485) (orb364694)

Description	Peptide to gp100 (476-485)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Val-Leu-Tyr-Arg-Tyr-Gly-Ser-Phe-Ser-Val
Formula	$C_{57}H_{83}N_{13}O_{15}$
Purity	> 95%
Source	Synthetic
MW	1190.38
Expiration Date	12 months from date of receipt.