

Product Datasheet

gp100 (570-579) (orb364695)

Description	Peptide to gp100 (570-579)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Leu-Ala-Asp-Thr-Asn-Ser-Leu-Ala-Val
Formula	$C_{41}H_{71}N_{11}O_{17}$
Purity	> 95%
Source	Synthetic
MW	990.09
Expiration Date	12 months from date of receipt.