

Product Datasheet

gp100 (639-647) (orb364698)

Description	Peptide to gp100 (639-647)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Leu-Pro-Arg-Ile-Phe-Cys-Ser-Cys
Formula	$C_{47}H_{79}N_{15}O_{11}S_2$
Purity	> 95%
Source	Synthetic
MW	1094.36
Expiration Date	12 months from date of receipt.