

Product Datasheet

PARP3 antibody (orb36484)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	This PARP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 451-481 amino acids of mouse PARP3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q8CFB8
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

