



Product Datasheet

S1PR3 antibody (orb36485)


Description	S1PR3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Sphingosine 1-phosphate receptor 3 protein (303-378AA)
Target	S1PR3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99500
Dilution Range	IHC-P:1:500-1000, IF:1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunofluorescent analysis of MCF-7 cell...



Immunohistochemistry of paraffin-embedde...