

## Product Datasheet

### GPR21 antibody (FITC) (orb36499)

|                        |                                                                              |
|------------------------|------------------------------------------------------------------------------|
| <b>Description</b>     | GPR21 antibody (FITC)                                                        |
| <b>Species/Host</b>    | Rabbit                                                                       |
| <b>Reactivity</b>      | Human                                                                        |
| <b>Conjugation</b>     | FITC                                                                         |
| <b>Immunogen</b>       | Recombinant Human Probable G-protein coupled receptor 21 protein (305-349AA) |
| <b>Target</b>          | GPR21                                                                        |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol               |
| <b>Form/Appearance</b> | Liquid                                                                       |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                |
| <b>Note</b>            | For research use only                                                        |
| <b>Isotype</b>         | IgG                                                                          |
| <b>Clonality</b>       | Polyclonal                                                                   |
| <b>Uniprot ID</b>      | <a href="#">Q99679</a>                                                       |
| <b>Dilution Range</b>  | ELISA                                                                        |
| <b>Expiration Date</b> | 12 months from date of receipt.                                              |