




Product Datasheet

GPC3 antibody (orb36547)

Description	Rabbit polyclonal antibody to GPC3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	GPC3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65563
Uniprot ID	P51654
NCBI	NP_004475.1 , NP_001158091.1 , NP_001158089.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



IHC-P of mouse stomach tissue lysates us...



Western Blot of HepG2 using Glypican 3 a...