




Product Datasheet

Metabotropic Glutamate Receptor 1 antibody (orb36550)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	GRM1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	132357
Uniprot ID	Q13255
NCBI	NP_001264994.1 , NP_001264996.1 , NP_001264995.1 , NP_001264993.1
Dilution Range	WB 1:50-1:100
Expiration Date	12 months from date of receipt.



WB analysis of human brain lysate using ...



Immunohistochemical staining of human st...