




Product Datasheet

IL1RA antibody (orb36621)

Description	Rabbit polyclonal antibody to IL1RA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	IL1RN
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20055
Uniprot ID	P18510
NCBI	NP_000568.1 , NP_776215.1 , NP_776213.1 , NP_776214.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50, IF/ICC: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of WiDr cell...



Immunohistochemical staining of paraffin...