

Product Datasheet

SCAPER antibody (Biotin) (orb36667)

Description	SCAPER antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human S phase cyclin A-associated protein in the endoplasmic reticulum protein (1109-1400AA)
Target	SCAPER
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BY12
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.