



Product Datasheet

HFE antibody (orb36682)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	HFE
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	40108
Uniprot ID	Q30201
NCBI	NP_620577.1 , NP_620580.1 , NP_000401.1 , NP_620573.1 , NP_620579.1 , NP_620572.1 , NP_620576.1 , NP_620578.1 , NP_620575.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of in Jurkat cell ...