



Product Datasheet

SPTLC3 antibody (orb36853)

Description	SPTLC3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Serine palmitoyltransferase 3 protein (1-58AA)
Target	SPTLC3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NUV7
Dilution Range	WB:1:1000-5000, IHC-P:1:200-500, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of HepG2 cell...



Western blot analysis of Mouse liver ti...