


Product Datasheet

FKBP14 antibody (orb36861)

Description	FKBP14 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP14 protein (20-211AA)
Target	FKBP14
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NWM8
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Mouse
stomach ...