



Product Datasheet

MCP1 antibody (orb36895)


Description	MCP1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human MCP1. Please contact us for the exact immunogen sequence. The peptide is available as orb374822.
Target	MCP1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	11 kDa
Uniprot ID	P13500
NCBI	002982, 002973, 31, 11
Entrez	6347
Dilution Range	IF/ICC: 1:50-500, IHC-P: 1:50-300, WB: 1:200-1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain...



IHC-P staining of rat brain tissue using...



IHC-P image of rat skin tissue using MCP...