




Product Datasheet

TIMP1 antibody (orb36909)

Description	Rabbit polyclonal antibody to TIMP1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	TIMP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23171
Uniprot ID	P01033
NCBI	NP_003245.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A2058 cel...



Immunohistochemical staining of paraffin...