

Product Datasheet

CCDC94 antibody (Biotin) (orb36944)

Description	CCDC94 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Coiled-coil domain-containing protein 94 protein (1-279AA)
Target	CCDC94
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BW85
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.