

Product Datasheet

NSUN3 antibody (orb36995)

Description	NSUN3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human tRNA (cytosine(34)-C(5))-methyltransferase, mitochondrial protein (1-200AA)
Target	NSUN3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H649
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.