




Product Datasheet

NFkB p65 antibody (orb37069)

Description	Rabbit polyclonal antibody to NFkB p65
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	RELA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60219
Uniprot ID	Q04206
NCBI	NP_001230913.1 , NP_068810.3 , NP_001138610.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of RELA antibody i...



Western blot analysis of RELA antibody i...