

Product Datasheet

Endothelin 1 antibody (orb37107)

Description	Rabbit polyclonal antibody to Endothelin 1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ECE1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	87164
Uniprot ID	P42892
NCBI	NP_001106819.1 , NP_001106818.1 , NP_001106820.1 , NP_001388.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

