



Product Datasheet

CD184 antibody (orb371619)

Description	Mouse monoclonal antibody to CD184
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P
Immunogen	CP-MAC-infected Sup-T1 cells
Target	CD184
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-4 µg/ml.
Isotype	Mouse IgG2a kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P61073
Entrez	7852
Dilution Range	Flow cytometry: Recommended dilution: 2-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of
human
CD184
positive
lymph...



Flow
cytometry
surface
staining
pattern ...