



Product Datasheet

Integrin alpha 4 antibody (orb371630)

Description	Mouse monoclonal antibody to Integrin alpha 4
Reactivity	Bovine, Canine, Equine, Feline, Human, Primate, Sheep
Conjugation	Unconjugated
Tested Applications	FA, FC, IHC-Fr
Target	CD49d
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P13612
Entrez	3676
Dilution Range	Functional application: In vitro T cell costimulation. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human CD49d positive lymph...



Flow cytometry surface staining pattern ...