

Product Datasheet

Keyhole limpet hemocyanin antibody (orb371639)

Description	Rabbit Polyclonal antibody to Keyhole limpet hemocyanin
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IL, WB
Immunogen	Purified keyhole limpet hemocyanin (KLH), whole molecule,
Form/Appearance	Liquid
Storage	Store at 4°C for 12-18 months. A preservative may be added for long time storage up to 2 years.
Note	For research use only
Application notes	Protein present in plant vascular tissue (xylem and vascular cambium) is detected by anti-KLH antibodies which might lead to false results in IL when using anti-peptide antibodies generated to KLH-conjugated peptide. Further information about it can be found here.
Clonality	Polyclonal
Purity	Immunogen affinity purified serum in PBS pH 7.4, conjugated to HRP.
MW	ca, 400 kDa/subunit
Uniprot ID	Q6KC56
Dilution Range	ELISA: 1:10000, WB: 1:10000
Expiration Date	12 months from date of receipt.