

Product Datasheet

GRO gamma/CXCL3 Antibody (orb371675)

Description

GRO gamma/CXCL3 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

E. coli-derived mouse GRO gamma recombinant protein (Position: S32-S100). Mouse GRO gamma shares 60.6% and 90.6% amino acid (aa) sequence identity with human and rat GRO gamma, respectively.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

10686 MW

Uniprot ID

[Q6W5C0](#)

Dilution Range

ELISA , 0.1-0.5µg/ml, Mouse, -Western blot, 0.1-0.5µg/ml, Mouse

Expiration Date

12 months from date of receipt.



WB analysis of GRO gamma using anti-GRO ...