


Product Datasheet

GIP Antibody (orb371680)

Description	GIP Antibody	100KD —
Species/Host	Rabbit	55KD —
Reactivity	Human, Rat	35KD —
Conjugation	Unconjugated	25KD —
Tested Applications	IHC, WB	15KD —
Immunogen	E. coli-derived human GIP recombinant protein (Position: Y52-Q93). Human GIP shares 92.9% and 95.2% amino acid (aa) sequence identity with mouse and rat GIP, respectively.	WB analysis of GIP using anti-GIP antibo...
Form/Appearance	Lyophilized	 IHC analysis of GIP using anti-GIP antib...
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	17108 MW	
Uniprot ID	P09681	
Dilution Range	Western blot, 0.1-0.5µg/ml, Human Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Human, Rat, By Heat	
Expiration Date	12 months from date of receipt.	