



Product Datasheet

CIITA Antibody (orb371730)

Description	CIITA Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	WB analysis of CIITA using anti-CIITA an...
Conjugation	Unconjugated	
Tested Applications	IHC, WB	IHC analysis of CIITA using anti-CIITA a...
Immunogen	E.coli-derived human CIITA recombinant protein (Position: E945-R1130). Human CIITA shares 84.4% amino acid (aa) sequence identity with mouse CIITA.	
Form/Appearance	Lyophilized	IHC analysis of CIITA using anti-CIITA a...
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	123514 MW	
Uniprot ID	P33076	
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Mouse, Rat, By Heat Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat	
Expiration Date	12 months from date of receipt.	