

Product Datasheet

DDAH1 Antibody (orb371734)


Description

DDAH1 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human DDAH1 (195-226aa QKALKIMQQMSDHRYSKLTVPDDIAANCIYLN), different from the related mouse and rat sequences by one amino acid.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31122 MW
Uniprot ID	O94760
Dilution Range	Western blot, 0.1-0.5µg/ml, Human, Mouse, RatImmunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Mouse, Rat, By Heat Immunocytochemistry/Immunofluorescence, 2µg/ml, Human Flow Cytometry, 1-3µg/1x10 ⁶ cells, Human
Expiration Date	12 months from date of receipt.



Flow
Cytometry
analysis of
A431 cells us...



WB analysis of
DDAH1 using
anti-DDAH1
an...



IF analysis of
DDAH1 using
anti-DDAH1
an...