



Product Datasheet

PF4 Antibody (orb371746)

Description	PF4 Antibody	100KD —
Species/Host	Rabbit	55KD —
Reactivity	Mouse	35KD —
Conjugation	Unconjugated	25KD — —
Tested Applications	ELISA, IHC, WB	15KD —
Immunogen	E. coli-derived mouse PF4 recombinant protein (Position: V30-S105). Mouse PF4 shares 74.3% and 86.8% amino acid (aa) sequence identity with human and rat PF4, respectively.	WB analysis of PF4 using anti-PF4 antibo...
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	IHC analysis of PF4 using anti-PF4 antib...
Clonality	Polyclonal	
MW	11243 MW	
Uniprot ID	Q9Z126	
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Mouse, Rat, By Heat ELISA , 0.1-0.5µg/ml, Mouse, - Western blot, 0.1-0.5µg/ml, Mouse	
Expiration Date	12 months from date of receipt.	