

## Product Datasheet

### Mouse BTN1A1 protein (orb371832)

<b>Description</b>	Recombinant mouse Butyrophilin Subfamily 1 Member A1 protein
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mMNaCl, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Application notes</b>	(C-6His)
<b>Protein Sequence</b>	Ala27-Trp247
<b>Source</b>	Human Cells
<b>Uniprot ID</b>	<b>Q62556</b>
<b>Expiration Date</b>	6 months from date of receipt.