

Product Datasheet

Mouse MCP1 protein (orb371839)

Description	Recombinant mouse C-C motif chemokine 2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris, 500mM NaCl, 10% glycerol, pH 7.4 .
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris, 500mM NaCl, 10% Glycerol, pH 7.4 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	(C-6His)
Protein Sequence	Gln24-Asn148
Source	Human Cells
Uniprot ID	P10148
Expiration Date	6 months from date of receipt.