

Product Datasheet

Mouse IL-23 protein (orb371847)

Description	Recombinant mouse Interleukin-23 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Val22-Ala196&Met23-Ser335
Source	Human Cells
Uniprot ID	Q9EQ14 , P43432
Expiration Date	6 months from date of receipt.