

Product Datasheet

Human ARG1 protein (orb371864)

Description	Recombinant human Arginase-1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 10mM Tris-HCl, 1mM 2-mercaptoethanol, 2mM MnCl ₂ , 150mM NaCl, pH 7.5.
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 1mM DTT, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	(C-6His)
Protein Sequence	Met1-Lys322
Source	Human Cells
Uniprot ID	P05089
Expiration Date	6 months from date of receipt.