

Product Datasheet

Mouse CCL3 protein (orb371883)

Description	Recombinant mouse C-C Motif Chemokine 3 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,5%Trehalose,1mM EDTA,pH8.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 5%Trehalose, 1mM EDTA, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Application notes	(N-6His)
Protein Sequence	Ala24-Ala92
Source	E.coli
Uniprot ID	P10855
Expiration Date	6 months from date of receipt.