

Product Datasheet

Mouse IL-27 protein (orb371900)

Description	Recombinant mouse Interleukin-27 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,2mM EDTA,5% Trehalose,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	(C-6His)
Protein Sequence	Tyr19-Pro228&Phe29-Ser234
Source	Human Cells
Uniprot ID	O35228, Q8K3I6
Expiration Date	6 months from date of receipt.